



The newest revision of the CO2 ejection system is now lighter, smaller and uses less pyrogen. We have made the e-match holders out of Delrin instead of stainless steel, so they can be used once. Of course, you can easily drill out the spent e-match and re-use the e-match holder.

The new CD4 also uses less pyrogen. This re-design produced better results vs the older CD3.

The new CD4 kits contain two sizes of CO2 cartridges- 16 gram and 38 gram. Over the last 15 years, these were the sizes most used by flyers. When calculating your ejection pressures needed, remember that each gram of CO2 produces 33 cubic inches of gas at seal level. Please review the sizing chart in the manual and on the website nickt written by John Rockdale.

SPECIFICATIONS

WEIGHT-	28.4 GRAMS	1.002 OUNCES
LENGTH FULLY ASSEMBLED-	2.52" LONG	
DIAMETER OF BODY-	.641"	
DIAMETER OF CAP- (MAX DIAMETER OF UNIT)	.808"	
CASING AND CAP-	6061-T6 ALUMINUM, ANODIZED RED.	
E-MATCH HOLDER-	DELIN, WHITE	
PLUNGER-	STAINLESS STEEL	